



SOP: LFA 00065

SOP For Backup/Restoring Analytical Data in E-Form

Alastair Sanderson

31/08/16

1 Objective

To be able to follow proper procedures to ensure accurate backup/restoration of analytical data from computers associated with the laboratory equipment.

2 Scope

Applicable to all analytical data from computers that are associated with laboratory equipment.

3 Responsibility

By: System Administrator

4 Accountability

Manager - Quality Control

5 Procedure

5.1 Analytical Record Archival

- 5.1. The data is acquired from the laboratory equipment through the computers running the software for processing and printing the chromatogram.
- 5.2. Analyst prints the processed data on paper.
- 5.3. Printed document will be examined and reviewed by the analyst and a reviewer.
- 5.4. The reviewed chromatogram with the original signature would serve as the master copy where all decisions pertaining to the material quality will be based.
- 5.5. The analytical record in paper would be archived and preserved as future reference for decision-making.

5.2 Data Back Up and Restore in Computer System

- 5.1. The computer system that is associated with the laboratory equipment will be a standalone system with limited storage space.
- 5.2. Due to the limited storage, only current data will be maintained in the computer system.
- 5.3. Periodic cleaning is necessary to delete old data from the computer and to ensure that there is available space for uninterrupted acquisition of data by the application software.
- 5.4. A maintenance log will be maintained by the system administrator to record the maintenance activity.
- 5.5. The system administrator will back up the current analytical data from the primary to the secondary hard disk in the computer system on a daily basis.
- 5.6. In the event that a system crash develops, the system administrator shall restore the system using the data from the secondary hard disk.
- 5.7. The restored data in the hard disk shall be maintained and can be produced on demand.

6 Abbreviations

SOP Standard Operating Procedure